

REALQUALITY RQ-HPV HR Multiplex — CE-IVD (Italy)

Detection of 14 high-risk HPV genotypes with specific 16/18 calls via Multiplex Real-Time PCR

Key Features

- Targets E6/E7 regions; minimizes false negatives linked to viral integration
- 14 HR coverage: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68
- Specific genotyping for HPV16 and HPV18
- Single-tube Multiplex RT-PCR; 5 µL extracted DNA input
- dUTP/UNG carry-over prevention; internal control and positive control included
- Harmonized thermal profile across REALQUALITY infectious/STI assays
- AB Genius Report for streamlined interpretation; GENEQUALITY® automation available
- Validated on common RT-PCR instruments (7500 Fast/Fast Dx, CFX96, LC480 II, RotorGene Q MDx)

Specimens

Extracted DNA from cervical/vaginal/urethral/oral/anal swabs; relevant biopsy/FFPE per IFU.

Ordering

- RQ974M — 50 tests
- RQ976M — 100 tests
- RQ97R4M — 50 tests (LightCycler® 480 II)
- RQ97R6M — 100 tests (LightCycler® 480 II)

For brochures/IFU and pricing, contact your local distributor.